

Cardio^{AI}

ECG annotation and interpretation software


Why Cardio.AI

It is widely accepted that long-term ECG diagnostics with the highest possible quality is hard to achieve, especially considering the current economical constraints. There are multiple factors that affect this quality, and the most prominent are:

- software that is used to clean up the record;
- technician's training level;
- time constraints.



Altogether these factors serve as the downward multipliers, reducing the chances to spot this single episode we might be looking for, affecting the efficiency of early detection and prevention. Multiple sources suggests the accuracy of the manual long-term ECG processing stays as low as 70%-75%.*



Pic.1. Only three episodes of AV-DISSOCIATION were detected which might be tricky to spot by using plain conventional Holter software.

* <https://jamanetwork.com/journals/jamainternalmedicine/article-abstract/2771093>

How we solve the problem

Cardio.AI solves this problem by introducing the AI pre-processing pipeline as the first step before the technician even attends the record. This mostly eliminates the human factor, improving the ability to spot arrhythmia or conduction issues. This particular AI separates sinus beats from everything else with very high accuracy, more than 95% in presence of the noise.

Technician's work transforms from manual processing (which is error-prone, see Pic.2) into validation processing in a much more relaxed setting as the AI already proposed interpretation (Pic.3). Technician needs only to validate and occasionally reclassify the strips.



Pic.2. Unannotated signal.



Pic.3. Automated annotation by Cardio.AI.

AI-assisted processing

That way we propose the best of both worlds - an entirely automatic pipeline but under full human control, allowing for occasional human intervention. AI's classification ability serves as the quality bottom line, supporting diagnostic companies with high-quality output.



Pic.4. Seeing is believing – all the findings at your fingertips.

To make things even more interesting, Cardio.AI already understands almost any possible arrhythmia or conduction issue, and the team constantly working towards reaching the full coverage, to include even rare arrhythmias, essentially a superhuman state on the ECG interpretation. The process of improving the AI is conceptually very simple, the team just adds the strips on which AI performs poorly to the training pipeline, eventually converging the quality to be above 90% (compared to human's average 74%) on any particular arrhythmia class, with some classes, like Atrial Fibrillation or Ventricular ectopic rhythms, to be above 98%.

Reporting

For physicians and attending clinicians all the above transforms into high-quality reports that frequently contains the references to the strips that would be otherwise hard or even impossible to spot with conventional Holter processing pipelines. This includes things like singular Type I AV-blocks, rhythm migration, nodal rhythm, etc.

The software resides completely in the cloud and is accessible from a number of supported web-browsers. This leads to another strong point – the ability for doctors to access the evidence online, from anywhere. Report becomes clickable for all the statistics and strips, leading to the signal in context. This allows for any exploration depth, a game-changer in some cases, particularly for patients with multiple arrhythmia types. It opens the internals of the record to the attending clinician, helping to build the trust between the doctor and a diagnostic center.

Public Demo

3, Republikas square, office 307, Riga, LV-1010, Latvia

Tel: +371 67 305 084 getintouch@xoresearch.com

Important

Date of Birth

1 Jan 1970
(50 yrs)

Gender

Unknown
gender

Uploaded by

Cardio.AI Support

Technician

Cardio.AI Support

Examination Time

26.11.2020 12:05 to
28.11.2020 12:07
2d 1m

Lead Configuration

I, II Leads

Heart Rate (Sinus)

Total Beats

221570

BPM

Max

147

18:20

Avg

79

Min

44

4:40

Table of Intervals

Max (ms)

Avg (ms)

Min (ms)

PQi

180

161

138

QT

410

361

321

QTc

417

408

400

ST - Segment

Max (mV)

Min (mV)

Dur

T-wave

I

0.05

0.01

4s

Λ

II

0.07

0.03

9s

Λ

Severity

Problem Area

Severity

Description

GENERAL

Automatism

✓

Low

Sinus Bradycardia, Sinus Tachycardia

Excitability

✓

Low

Wandering Sinus Pacemaker Within The Sinus Node, Atrial Ectopic Rhythm, Upper Atrial Rhythm

Conductivity

✓

Important

First Degree AV Block, Second Degree AV Block Type I, AV Dissociation With Interference

Pre-excitation

✓

Ok

ECTOPIA

Supraventricular

✓

Low

Atrial Escape Beat

Ventricular

✓

Ok

Tachycardia

✓

Ok

Fibrillation/Flutter

✓

Ok

STEMI

ST depression/elevation

✓

Ok

T-wave

✓

Ok

PQRST

PQ interval/segment

✓

Ok

QRS complex

✓

Ok

QTc interval

✓

Ok

Pic.5. Decision-supporting reports.