Process Development

Services


Equipment with guaranteed performances

- Feasibility studies and process line definition
- Process modeling, process optimization, robustness
- Process validation and scale-up
- Sample production (on customer site or at our premises)

A Unique Value Proposition

Each year, Novasep's team develops about 30 processes for the Industrial Biotechnology industries. In particular, our team can:

- · Define the most appropriate process line
- Determine your scale-up parameters
- Simulate your process performances
- Produce your samples from kg to ton scale
- Assure clean and safe procedures
- Guarantee the performances of your process



Bio-Industries Bio-based chemicals Amino acids Organic acids Antibiotics Vitamins

0000





Oligosaccharides Sweeteners Fatty Acids Antioxidants Proteins



Novasep Enables your Targeted Product Development

Our Facilities

- 50 mobile pilot units
- · A global presence with locations in France, USA and China
- · Possible piloting or pilot line building on any customer site worldwide
- · A team of 25 dedicated and highly skilled people








400 m² pilot hall and Labs, China

Saint-Maurice-De-Beynost, France

Our technologies



Analytical Development Capabilities

Novasep applies the following analytical tools:

- · HPLC/UPLC with UV-Vis, DAD, ELSD and RI detectors
- · Ion chromatography (Conductimetry detector)
- Gas chromatography and Headspace GC
- Spectrophotometry
- SDS-PAGE electrophoresis
- Mass spectrometry
- ELISA development




Our Key Technologies Available at Pilot Scale

Batch Chromatography and Ion Exchange

- Columns range from 20 L to 800 L
- Fully automated skids available

Continuous Chromatography and Ion Exchange

- SSMB (Sequential Simulated Moving Bed) skids available with 2 to 12 columns
- Continuous Ion Exchange skids with up to 10 columns







Simulation software

SSMB Skid

Membrane Filtration

A wide range of membranes can be used:

- Organic membranes for Nanofiltration and Reverse osmosis from 0.3 m² to 240 m²
- Mineral membranes for Micro and Ultrafiltration from 0.2 m² to 10 m²

Novasep uses its proprietary ceramic membranes Kerasep®





Ultrafiltration Skid and Kerasep® module


Electrodialysis (ED)

Membranes and systems for:

- Conventional ED
- Bipolar Membrane ED from 2 to 4 m²




Electrodialysis Skid





Evaporation Skid

Case studies: process lines fully developed by our teams



Some examples of our sample production capacities:



· Batch chromatography: 200 g of high purity plant extract

- Electrodialysis: 100 kg of demineralized lactose (D50, D70 & D90)
- · SSMB continuous chromatography: 500 kg of purified oligosaccharides
- Filtration: up to 6 m3/day for broth filtration



EUROPE - Novasep Phone: +33 4 72 01 27 27 industrialbiotech@novasep.com

INDIA - Novasep India Phone: +91 900 84 77 311 industrialbiotech.india@novasep.com NORTH AMERICA - Novasep Inc Phone: +1 610 494 0447 industrialbiotech.north-am@novasep.com

RUSSIA - Novasep Russia Phone: +7 (8652) 316 645 industrialbiotech.russia@novasep.com ASIA - Novasep Asia Phone: +86 21 6045 1600 industrialbiotech.asia@novasep.com

REST OF THE WORLD - Novasep Phone: +33 4 72 01 27 27 industrialbiotech@novasep.com

Services and technologies for the life science industries

www.novasep.com novasep@novasep.com

JAPAN - Novasep Japan K.K. Phone: +81 (0)3 3221 1960 japan@novasep.com

© 2014 Groupe Novasep SAS. Groupe Novasep SAS and/or its affiliates (hereafter Groupe Novasep) own or license the copyrights, trademarks, names, logos, and other rights to the information on this brochure. No right or license is granted to any other party. Any unauthorized use without the express prior written consent of Groupe Novasep is prohibited. Disclaimers: The information contained in this brochure: Are provided "as is", for informational purposes only without any representation or warranty of accuracy or completeness. In no event will Groupe Novasep be liable to any party for any damages or any other losses from access or reliance upon any information contained in this brochure. 033-001B