

# **Single Patient Use Transverse Cutter and Reload**

#### **Major Features**

#### **Great Operation Experience**

- O Curved small head, easy to enter deep pelvic under direct vision
- O Intermediate position, easy to adjust cutting and stapling location
- o Separate closing and firing handles, clear operation processes

#### **Great Efficiency**

- o Short reload width, long cutting line
- O Simultaneous cutting and stapling, simplified procedures
- O Quick and convenient operation

## **Great Reliability**

- o Stapling line 4.5mm longer than cutting line
- Tactile and audible feedback in operation process
- o Positioning pin, ensured stapling of all target tissue
- o Parallel closure, identical formation of all staples

Golden Stapler Surgical Co., Ltd.









**Advancing Surgical Technology** to Benefit More People









#### **Curved Small Head, 62mm Width**

Easier for vertical placement into deep pelvic Procedure under surgeon's direct vision 4 rows of staples released during firing, tissue in the middle to be cut Stapling while cutting, easy and convenient Reduce surgery processes, more effective



#### **Positioning Pin**

Close and lock target tissue, prevent tissue slippage Assist parallel closure Uniform B-shape formation of all staples Great stapling and hemostasis effect



#### **Curved Reload**

Short in width, long in curved cutting line Compress tissue up to 40mm Stapling line 4.5mm longer than cutting line Great stapling and hemostasis in cutting line's end



### **Jaw Opening 19mm Wide**

Big opening, easy to put tissue inside





### **Separate Closing and Firing Handles**

Independent closing and firing processes Prevent miss-operation



## **Intermediate Position**

Easy to move the anvil to deep pelvic Easier to adjust cutting and stapling location



#### **Cutting Washer**

Tactile and audible feedback during firing Clear suggestion of operation process

Code	Color	Cutting Line	Open Staple Height	Staple Rows	Closed Staple Height	Staple Qty	Tissue Thickness
CCSA40G	Green	40mm	4.5mm	4	2.0mm	46	Thick

<sup>\*</sup> Staple from titanium material

Golden Stapler Surgical Co., Ltd.

Add: Unit B, 26 Huashan Zhong Road, Changzhou, Jiangsu, China, 213022 Tel: +86 519 85266817 Fax: +86 519 85266821

Web: http://www.goldenstapler.com