

Single Patient Use Linear Stapler and Reload

Available of Yellow Reload and 45mm Hand Piece

- White, blue, yellow & green reloads to accommodate to various tissue thicknesses
- 45mm linear stapler, ideal for double stapling technique in low rectal sphincter-preserving procedure

Social States States 45m

Major Features

- Stapling and closure of intestinal, pulmonary, gastric and vascular incision
- Available of 45 and 60mm models with blue, yellow and green reloads;
 30mm model with blue and green reloads
- Unique 30mm vessel stapler with white reload of 3 staple rows, specialized for the best vessel stapling
- One single hand piece can be fired 8 times

Golden Stapler Surgical Co., Ltd.









Advancing Surgical Technology to Benefit More People



Single Patient Use Linear Stapler and Reload

Positioning Pin

- Avoid slippage of clamped tissue
- Assist parallel closure
- Ensure complete suturing of the target tissue

Parallel Closure

- Even compression on target tissue
- · Provide uniform formation height of all staples
- · Effective suturing & homeostasis



- White reload, 3 rows of staples, exclusive for vessel stapling
- Best stapling and hemostasis effect

Blue, Yellow, Green & White Reloads

- Accommodate to various tissue thicknesses
- · Effective stapling and hemostasis

3 Different Sizes 30, 45

- · Accomodate various GI tract

& 60mm

- · Economic and easy to use

Intermediate Position

- · Adjust suturing position on target tissue
- · Ensure suturing on target position

Separate Closure & **Firing Handles**

- · Clear surgical procedures
- · Safe & easy to use





Tactile & Audible Responses

 Safety mechanism for empty reload to prevent misfiring

| Code | Colour | Staple Rows | Staple Qty | Closed Staple Height (MM) | Stapling Length (MM) | Applicable Tissue Thickness | Reload Code |
|--------|--------|-------------|------------|---------------------------------|-------------------------|--------------------------------|-------------|
| LSC30W | White | 3 | 23 | 1.0 | 30 | Blood vessel, thin tissue | LSCR30W |
| LSC30B | Blue | 2 | 11 | 1.5 | 30 | Normal thick tissue | LSCR30B |
| LSC30G | Green | 2 | 11 | 2.0 | 30 | Thick tissue | LSCR30G |
| LSC45B | Blue | 2 | 15 | 1.5 | 45 | Normal thick tissue | LSCR45B |
| LSC45Y | Yellow | 2 | 15 | 1.8 | 45 | Medium thick tissue | LSCR45Y |
| LSC45G | Green | 2 | 15 | 2.0 | 45 | Thick tissue | LSCR45G |
| LSC60B | Blue | 2 | 21 | 1.5 | 60 | Normal thick tissue | LSCR60B |
| LSC60Y | Yellow | 2 | 21 | 1.8 | 60 | Medium thick tissue | LSCR60Y |
| LSC60G | Green | 2 | 21 | 2.0 | 60 | Thick tissue | LSCR60G |

^{*} Staple from titanium material

Golden Stapler Surgical Co., Ltd.

Add: Building 7A, Jiangsu Wujin Sci-Tech Pioneer Park, 256 Mid Mingxin Road, Hutang Town, Wujin District, Changzhou, Jiangsu Province, China 213164

Tel: +86 519 85070251 Fax: +86 519 85070252

Web: http://www.goldenstapler.com